

CUTIMED SILTEC PLUS FOAM DRESSING



MAIN SPECIFICATIONS

- Cutimed Siltec PLUS through enhanced grip, remains in position while applying the secondary fixation, ensuring atraumatic and painless removal.
- Superabsorbent polyurethane foam dressing, with a silicone contact layer and soft grip, without adhesive border.
- Moderate to high exudate levels.
- 5 mm thickness.
- Padded protective effect.
- Siltec PLUS also available for injured heel.

| Code | Description | Selling unit |
|----------|-----------------------------------|--------------|
| 72683-05 | CUTIMED SILTEC PLUS HEEL 16x24 CM | CX.5 |
| 72683-04 | CUTIMED SILTEC PLUS 20x20 CM | CX.5 |
| 72683-03 | CUTIMED SILTEC PLUS 15x15 CM | CX.10 |
| 72683-02 | CUTIMED SILTEC PLUS 10x20 CM | CX. 10 |
| 72683-01 | CUTIMED SILTEC PLUS 10x10 CM | CX. 10 |
| 72683-00 | CUTIMED SILTEC PLUS 5x6 CM | CX. 10 |

FACTORMED
Produtos Médicos & Hospitalares