

CUTIMED SILTEC L (LIGHT) FOAM DRESSING



MAIN SPECIFICATIONS

- Superabsorbent polyurethane foam dressing, with silicone contact layer, without adhesive edge.
- Low to moderate exudate levels.
- High comfort and safe management of exudate, even under compression.

Code	Description	Selling unit
72630-02	CUTIMED SILTEC L 15x15 CM	CX. 10
72630-01	CUTIMED SILTEC L 10x10 CM	CX. 10
72630-00	CUTIMED SILTEC L 5x6 CM	CX. 10