

CUTIMED SILTEC HEEL 2D FOAM DRESSING



MAIN SPECIFICATIONS

- Superabsorbent polyurethane foam dressing with silicone contact layer for the prevention and treatment of heel wounds.
- Moderate to high exudate levels.
- Two adhesion levels.
- Additional cushioned protective effect (7 mm), helps to reduce the risks from forces of pressure, friction and shear.

Code	Description	Selling unit
1254006	CUTIMED SILTEC HEEL 2D 16x24 CM	CX.5