

CUTIMED SILTEC SACRUM FOAM DRESSING



MAIN SPECIFICATIONS

- Superabsorbent polyurethane foam dressing with silicone contact layer and silicone adhesive edge for optimum adaptation to the sacred area.
- Moderate to high exudate levels.
- Anatomically designed for perfect adaptation to the perianal area.
- Effective exudate management, even in complex situations such as pressure ulcers in the sacred area.
- Trip 3 cm around the perimeter for a better fixation and greater safety of treatment, in particular for incontinent patients.

Code	Description	Selling unit
72647-01	CUTIMED SILTEC SACRUM 23x23 CM	CX.5
1254005	CUTIMED SILTEC SACRUM 17,5x17,5 CM	CX.5