

# CUTIMED ALGINATE CALCIUM ALGINATE DRESSING



## MAIN SPECIFICATIONS

- The Cutimed Alginate is a calcium sodium alginate dressing (80% calcium and 20% sodium) for the management of highly exuding wounds.
- High absorption capacity.
- Quick gelification.
- It absorbs exudate to the gel, reducing the risk of leakage and maceration.
- The soft, cohesive gel creates and maintains a moist wound environment.
- Good wet integrity: allows quick and easy removal as a single piece.
- Soft and malleable: perfect adaptation to the morphology of the wound.

Code	Description	Selling unit
72634-03	CUTIMED ALGINATE 2,5x30 CM	CX.5
72634-02	CUTIMED ALGINATE 10x20 CM	CX. 10
1254016	CUTIMED ALGINATE 10x10 CM	CX. 10
72634-00	CUTIMED ALGINATE 5x5 CM	CX. 10



[www.factormed.pt](http://www.factormed.pt)