

CUTIMED ALGINATE CALCIUM ALGINATE DRESSING



MAIN SPECIFICATIONS

- The Cutimed Alginate is a calcium sodium alginate dressing (80% calcium and 20% sodium) for the management of highly exuding wounds.
- High absorption capacity.
- Quick gelification.
- It absorbs exudate to the gel, reducing the risk of leakage and maceration.
- The soft, cohesive gel creates and maintains a moist wound environment.
- Good wet integrity: allows quick and easy removal as a single piece.
- Soft and malleable: perfect adaptation to the morphology of the wound.

| Code | Description | Selling unit |
|----------|----------------------------|--------------|
| 72634-03 | CUTIMED ALGINATE 2,5x30 CM | CX.5 |
| 72634-02 | CUTIMED ALGINATE 10x20 CM | CX. 10 |
| 1254016 | CUTIMED ALGINATE 10x10 CM | CX. 10 |
| 72634-00 | CUTIMED ALGINATE 5x5 CM | CX. 10 |

FACTORMED
Produtos Médicos & Hospitalares