

# STERILIZED SHEET FOR BURNED



## MAIN SPECIFICATIONS

- Sterile sheet for burns, in nonwoven hydrophobic of low adhesion, 100% polypropylene.
- To involve the wounded burned from the scene of the accident to the hospital.
- Large dimensions to fully engage the patient.
- Blue color.
- Weight 40g / m<sup>2</sup>.

Code	Description	Selling unit
39016025	STERILIZED SHEET FOR BURNED 150X240 CM	CX. 1
39016032	STERILIZED SHEET FOR BURNED 150X180 CM	CX. 1