

COLOHESIVE PLUS COLOSTOMY POUCH 1 PIECE CLOSED OPAQUE



MAIN SPECIFICATIONS

- Great convenience and discretion.
- High skin protection.
- With Stomahesive protection: In contact with the peristomal skin, offers maximum protection while having a healing effect in cases of skin abrasion.
- Incorporates impregnated charcoal filter with good deodorizing ability, ideal for device discretion.
- Covered with a lacy/porous texture to allow the absorption of perspiration, avoiding the unpleasant effect of "paste" the abdomen.
- Opaque.

| Code | Description | Selling unit |
|-------|---|--------------|
| 62450 | COLOHESIVE PLUS 1 PIECE CLOSED OPAQUE 45 MM | CX.30 |
| 62449 | COLOHESIVE PLUS 1 PIECE CLOSED OPAQUE 38 MM | CX. 30 |
| 62448 | COLOHESIVE PLUS 1 PIECE CLOSED OPAQUE 32 MM | CX.30 |
| 62451 | COLOHESIVE PLUS 1 PIECE CLOSED OPAQUE 50 MM | CX. 30 |

