

COLOHESIVE PLUS COLOSTOMY POUCH 1 PIECE CLOSED OPAQUE



MAIN SPECIFICATIONS

- Great convenience and discretion.
- High skin protection.
- With Stomahesive protection: In contact with the peristomal skin, offers maximum protection while having a healing effect in cases of skin abrasion.
- Incorporates impregnated charcoal filter with good deodorizing ability, ideal for device discretion.
- Covered with a lacy/porous texture to allow the absorption of perspiration, avoiding the unpleasant effect of "paste" the abdomen.
- Opaque.

Code	Description	Selling unit
62450	COLOHESIVE PLUS 1 PIECE CLOSED OPAQUE 45 MM	CX.30
62449	COLOHESIVE PLUS 1 PIECE CLOSED OPAQUE 38 MM	CX. 30
62448	COLOHESIVE PLUS 1 PIECE CLOSED OPAQUE 32 MM	CX.30
62451	COLOHESIVE PLUS 1 PIECE CLOSED OPAQUE 50 MM	CX. 30

