

# ESTEEM SYNERGY COLOSTOMY POUCH 2 PIECES STANDARD CLOSED OPAQUE



## MAIN SPECIFICATIONS

- The Esteem Synergy ostomy system features an innovative Adhesive Coupling Technology. Instead of plastic rings, a foam adhesive allows the pouch to be securely attached to the skin barrier on a clear plastic “Landing Zone” flange. The result is a low profile, flexible, two-piece system.
- Pouch is designed for single use.
- Simply discard when full.

Code	Description	Selling unit
403946	ESTEEM SYNERGY 2 PIECES STANDARD CLOSED OPAQUE 61 MM	CX.30
405437	ESTEEM SYNERGY 2 PIECES STANDARD CLOSED OPAQUE 48 MM	CX.30
405436	ESTEEM SYNERGY 2 PIECES STANDARD CLOSED OPAQUE 35 MM	CX.30