

ESTEEM SYNERGY COLOSTOMY POUCH 2 PIECES STANDARD CLOSED TRANSPARENT



MAIN SPECIFICATIONS

- The Esteem Synergy ostomy system features an innovative Adhesive Coupling Technology. Instead of plastic rings, a foam adhesive allows the pouch to be securely attached to the skin barrier on a clear plastic “Landing Zone” flange. The result is a low profile, flexible, two-piece system.
- Pouch is designed for single use.
- Simply discard when full.

Code	Description	Selling unit
405439	ESTEEM SYNERGY 2 PIECES STANDARD CLOSED TRANSPARENT 48 MM	CX. 30
405438	ESTEEM SYNERGY 2 PIECES STANDARD CLOSED TRANSPARENT 35 MM	CX.30
405440	ESTEEM SYNERGY 2 PIECES STANDARD CLOSED TRANSPARENT 61 MM	CX. 30