

CONSECURA AUTOLOCK COLOSTOMY POUCH 2 PIECES CLOSED TRANSPARENT WITH FILTER



MAIN SPECIFICATIONS

- The ostomy system Consecura Autolock is an automatic fastening technology.
- Soft touch, the coupling mechanism between the pouch and the protective is fixed safely, minimizing the pressure in the abdomen.
- Easy to use and is very useful for people with poor coordination.
- An audible "click" of the engagement ring, confirms secure connection to the skin protector.
- The coupling allows a 360° bag rotation, remaining fixed in place, which offers convenience and ease of product use.
- Comfort panels, smooth and double size are designed to be quiet and odor-proof, helping discretion maintenance.

Code	Description	Selling unit
401313FI	CONSECURA AUTOLOCK 2 PIECES CLOSED TRANSPARENT WITH FILTER 70 MM	CX. 30

401311FI	CONSECURA AUTOLOCK 2 PIECES CLOSED TRANSPARENT WITH FILTER 45 MM	CX. 30
----------	---	--------

401312FI	CONSECURA AUTOLOCK 2 PIECES CLOSED TRANSPARENT WITH FILTER 57 MM	CX. 30
----------	---	--------
