

COMBIHESIVE IIS COLOSTOMY POUCH 2 PIECES CLOSED WITH FILTER



MAIN SPECIFICATIONS

- Device with two pieces connected, comfortable, flexible and secure.
- This bag consists of two separate elements: a support and a pouch. The bag can be changed without removing the support. The bag is fastened to the support through interlocking rings. we can hear a click when the bag is securely attached to the support.
- It offers a full range of moldable bags to best fit to the stoma, without having to cut the bag, allowing for ideal sealing equipment.
- The bag may be easily changed without removing the support.
- It provides safety, comfort and skin protection to the stoma.
- Filter for better control of gases and smells.
- Texture more comfortable to the touch and more discreet.
- Opaque bag with anatomical shape for better distribution of the content and ease of use.
- Internal structure with eight retaining teeth for secure coupling bag/dressing. Tactile.

Code

Description

Selling unit

402525	COMBIHESIVE IIS 2 PIECES CLOSED WITH FILTER 70 MM	CX.30
402524	COMBIHESIVE IIS 2 PIECES CLOSED WITH FILTER 57 MM	CX.30
402523	COMBIHESIVE IIS 2 PIECES CLOSED WITH FILTER 45 MM	CX.30
402522	COMBIHESIVE IIS 2 PIECES CLOSED WITH FILTER 38 MM	CX.30
402521	COMBIHESIVE IIS 2 PIECES CLOSED WITH FILTER 32 MM	CX.30