

# ESTEEM SYNERGY ILEOSTOMY 1 PIECE STANDARD DRAINABLE TRANSPARENT



## MAIN SPECIFICATIONS

- The Esteem Synergy ostomy system features an innovative Adhesive Coupling Technology. Instead of plastic rings, a foam adhesive allows the pouch to be securely attached to the skin barrier on a clear plastic “Landing Zone” flange. The result is a low profile and flexible system.
- Pouch is designed for single use.
- Simply discard when full.

Code	Description	Selling unit
405406	ESTEEM SYNERGY 1 PIECE STANDARD DRAINABLE TRANSPARENT 89 MM	CX. 10 + 2 CLIPS
405405	ESTEEM SYNERGY 1 PIECE STANDARD DRAINABLE TRANSPARENT 61 MM	CX. 10 + 2 CLIPS
405404	ESTEEM SYNERGY 1 PIECE STANDARD DRAINABLE TRANSPARENT 48 MM	CX. 10 + 2 CLIPS