

# CONSECURA AUTOLOCK ILEOSTOMY POUCH 1 PIECE DRAINABLE TRANSPARENT



## MAIN SPECIFICATIONS

- The ostomy system Consecura Autolock is an automatic fastening technology.
- Soft touch, the coupling mechanism between the pouch and the protective is fixed safely, minimizing the pressure in the abdomen.
- Easy to use and is very useful for people with poor coordination.
- An audible "click" of the engagement ring, confirms secure connection to the skin protector.
- The coupling allows a 360° bag rotation, remaining fixed in place, which offers convenience and ease of product use.
- Comfort panels, smooth and double size are designed to be quiet and odor-proof, helping discretion maintenance.

Code	Description	Selling unit
401443FI	CONSECURA AUTOLOCK 1 PIECE DRAINABLE TRANSPARENT 100 MM	CX. 10 + 1 CLIP

401327FI	CONSECURA AUTOLOCK 1 PIECE DRAINABLE TRANSPARENT 45 MM	CX. 10 + 1 CLIP
----------	---	-----------------

---

401441FI	CONSECURA AUTOLOCK 1 PIECE DRAINABLE TRANSPARENT 57 MM	CX. 10 + 1 CLIP
----------	---	-----------------

---