

# COMBIHESIVE IIS ILEOSTOMY POUCH 2 PIECES DRAINABLE TRANSPARENT



## MAIN SPECIFICATIONS

- Device with two pieces connected, comfortable, flexible and secure.
- This bag consists of two separate elements: a support and a pouch. The bag can be changed without removing the support. The bag is fastened to the support through interlocking rings. we can hear a click when the bag is securely attached to the support.
- It offers a full range of moldable bags to best fit to the stoma, without having to cut the bag, allowing for ideal sealing equipment.
- The bag may be easily changed without removing the support.
- It provides safety, comfort and skin protection to the stoma.
- Texture more comfortable to the touch and more discreet.
- Transparent bag with anatomical shape for better distribution of the content and ease of use.
- Internal structure with eight retaining teeth for secure coupling bag/dressing. Tactile.

Code	Description	Selling unit
400970	COMBIHESIVE IIS 2 PIECES DRAINABLE TRANSPARENT 70 MM	CX. 10 + 2 CLIPS

400969	COMBIHESIVE IIS 2 PIECES DRAINABLE TRANSPARENT 57 MM	CX. 10 + 2 CLIPS
400968	COMBIHESIVE IIS 2 PIECES DRAINABLE TRANSPARENT 45 MM	CX. 10 + 2 CLIPS
400967	COMBIHESIVE IIS 2 PIECES DRAINABLE TRANSPARENT 38 MM	CX. 10 + 2 CLIPS
400971	COMBIHESIVE 2 PIECES DRAINABLE TRANSPARENT 100 MM	CX. 10 + 2 CLIPS