

CONSECURA AUTOLOCK ACCUSEAL UROSTOMY POUCH TRANSPARENT



MAIN SPECIFICATIONS

- The system closes and locks automatically with a single light touch. An audible click confirms the attachment of skin barrier to pouch. Pouches can be rotated or repositioned while locked in place. Tap allows for draining of urine.

Code	Description	Selling unit
401452FI	CONSECURA AUTOLOCK ACCUSEAL TRANSPARENT 35 MM	CX. 10 + 2 adapt.
401453FI	CONSECURA AUTOLOCK ACCUSEAL TRANSPARENT 45 MM	CX. 10 + 2 adapt.
401460FI	CONSECURA AUTOLOCK ACCUSEAL TRANSPARENT 57 MM	CX. 10 + 2 adapt.
401461FI	CONSECURA AUTOLOCK ACCUSEAL TRANSPARENT 70 MM	CX. 10 + 2 adapt.