

COMBIHESIVE IIS ACCUSEAL UROSTOMY POUCH



MAIN SPECIFICATIONS

- Comprehensive product line with the latest skin barrier technologies. The plastic coupling ring allows the pouch and skin barrier to snap together for audible security. Combihesive IIS urostomy pouches have an extended anti-reflux valve, comfort panels and a choice of two types of drainage taps.
- Accuseal drainage system.

Code	Description	Selling unit
400988	COMBIHESIVE IIS ACCUSEAL 32 MM	CX. 10 + 2 adapt.
400989	COMBIHESIVE IIS ACCUSEAL 38 MM	CX. 10 + 2 adapt.
400990	COMBIHESIVE IIS ACCUSEAL 45 MM	CX. 10 + 2 adapt.
400991	COMBIHESIVE IIS ACCUSEAL 57 MM	CX. 10 + 2 adapt.