

# COMBIHESIVE IIS UROSTOMY POUCH



## MAIN SPECIFICATIONS

- Comprehensive product line with the latest skin barrier technologies. The plastic coupling ring allows the pouch and skin barrier to snap together for audible security. Combihesive IIS urostomy pouches have an extended anti-reflux valve, comfort panels and a choice of two types of drainage taps.

Code	Description	Selling unit
402548	COMBIHESIVE IIS 32 MM	CX. 10 + 2 adapt.
402549	COMBIHESIVE IIS 38 MM	CX. 10 + 2 adapt.
402550	COMBIHESIVE IIS 45 MM	CX. 10 + 2 adapt.
402551	COMBIHESIVE IIS 57 MM	CX. 10 + 2 adapt.
402552	COMBIHESIVE IIS 70 MM	CX. 10 + 2 adapt.