

ESTEEM SYNERGY STOMAHESIVE DRESSING



MAIN SPECIFICATIONS

- Esteem synergy is an innovative device for stomas with adhesive fastening technology. Instead of plastic rings, a patented adhesive foam allows the adhesion of the bag in a transparent "association board". For this reason the 2-piece system is thin and flexible. The Stomahesive skin barriers are designed to be gentle on the peristomal skin.

Code	Description	Selling unit
405484	ESTEEM SYNERGY STOMAHESIVE DRESSING 13/35 MM	CX. 5
405485	ESTEEM SYNERGY STOMAHESIVE DRESSING 13/48 MM	CX. 5
405486	ESTEEM SYNERGY STOMAHESIVE DRESSING 13/61 MM	CX. 5
405487	ESTEEM SYNERGY STOMAHESIVE DRESSING 13/89 MM	CX. 5