


ESTEEM SYNERGY MOLDABLE FLEXIBLE DRESSING



MAIN SPECIFICATIONS

- Esteem synergy is an innovative device for stomas with adhesive fastening technology. Instead of plastic rings, a patented adhesive foam allows the adhesion of the bag in a transparent "association board". For this reason the 2-piece system is thin and flexible. The mouldable adhesive provides greater protection of the skin, in comparison with traditional plates.

Code	Description	Selling unit
409271	ESTEEM SYNERGY MOLDABLE FLEXIBLE DRESSING 13/22 MM POUCH 35 MM	CX. 5
409272	ESTEEM SYNERGY MOLDABLE FLEXIBLE DRESSING 22/33 MM POUCH 48 MM	CX. 5
409273	ESTEEM SYNERGY MOLDABLE FLEXIBLE DRESSING 33/45 MM POUCH 61 MM	CX. 5