

# CONSECURA AUTOLOCK FLEXIBLE DRESSING



## MAIN SPECIFICATIONS

- The ConSecura system features automatic closing with a simple operation. An audible "click" confirms the attachment of the pouch to the skin barrier. The skin barrier helps protecting the skin around the stoma.

Code	Description	Selling unit
401429FI	CONSECURA AUTOLOCK FLEXIBLE DRESSING 35 MM	CX. 5
401430FI	CONSECURA AUTOLOCK FLEXIBLE DRESSING 45 MM	CX. 5
401431FI	CONSECURA AUTOLOCK FLEXIBLE DRESSING 57 MM	CX. 5
401432FI	CONSECURA AUTOLOCK FLEXIBLE DRESSING 70 MM	CX. 5