


# COMBIHESIVE IIS DURAHESIVE CONVEX DRESSING



## MAIN SPECIFICATIONS

- Square plate with rounded, supple and soft corners, containing the Durahesive protective skin barrier, laminated onto a convex disk with a central hole and a low-profile flange. The surface of the board is laminated with a thin, clear polyethylene film. Upon the polyethylene film surface around the flange, the adhesive support is placed. The adhesive plate is protected by a silicone paper that must be removed at the time of use.

Code	Description	Selling unit
125028	COMBIHESIVE IIS DURAHESIVE CONVEX DRESSING 19/45 MM	CX. 5
125029	COMBIHESIVE IIS DURAHESIVE CONVEX DRESSING 22/45 MM	CX. 5
125030	COMBIHESIVE IIS DURAHESIVE CONVEX DRESSING 25/45 MM	CX.5
125031	COMBIHESIVE IIS DURAHESIVE CONVEX DRESSING 28/45 MM	CX. 5
125032	COMBIHESIVE IIS DURAHESIVE CONVEX DRESSING 32/45 MM	CX. 5
125033	COMBIHESIVE IIS DURAHESIVE CONVEX DRESSING 35/45 MM	CX. 5