


# COMBIHESIVE IIS LITTLE ONES PEDIATRIC DRESSING



## MAIN SPECIFICATIONS

- Stomahesive flexible children board, developed to protect the skin around the children stomas against irritation caused by effluent. Combihesive Iis Little Ones plate is composed of gelatin, pectin, sodium carboxymethylcellulose, polyisobutylene, polyethylene, polyester and acrylic adhesive.

Code	Description	Selling unit
125016	COMBIHESIVE IIS LITTLE ONES DRESSING 32 MM	CX. 5