

# LEUKOMED T PLUS DRESSING



## MAIN SPECIFICATIONS

- Postsurgical dressing with high absorption pad, which allows the patient to continue their daily activities. It is ideal for protecting wounds with moderate exudate levels, such as post-surgical wounds, lacerations, cuts, abrasions or minor burns.
- Breathable, reduces the risk of maceration. It increases the dressing durability.
- Waterproof, facilitates personal hygiene.
- Ultra-thin, soft, flexible film that adapts to body contours like a second skin without restricting movement.
- Sterile. Barrier against bacteria and external contaminants while reducing the risk of infection.
- Secure adhesion and high skin tolerance that holds the patch in place. Greater comfort for the patient.
- Highly absorbent pad with a low adhesion layer that contact does not interfere with the healing process and minimizes pain on removal. Keeps the wound dry.
- Rounded edges that prevent premature wear.

Code	Description	Selling unit
72382-00	LEUKOMED T PLUS 5 CM X 7,2 CM (PAD 2,8 CM X 4 CM)	CX. 50
72382-01	LEUKOMED T PLUS 8 CM X 10 CM (PAD 4 CM X 6 CM)	CX. 50
72382-02	LEUKOMED T PLUS 8 CM X 15 CM (PAD 4 CM X 11 CM)	CX. 50

72382-14	LEUKOMED T PLUS 10 CM X 20 CM (PAD 5 CM X 16 CM)	CX. 50
72382-03	LEUKOMED T PLUS 10 CM X 25 CM (PAD 5 CM X 20,5 CM)	CX. 50
72382-04	LEUKOMED T PLUS 10 CM x 30 CM (PAD 5 CM X 25 CM)	CX. 50
72382-05	LEUKOMED T PLUS 10 CM x 35 CM (PAD 5 CM X 30 CM)	CX. 50
72382-06	LEUKOMED T PLUS 5 CM X 7,2 CM (PAD 2,8 CM X 4 CM)	CX. 5
72382-07	LEUKOMED T PLUS 8 CM X 10 CM (PAD 4 CM X 6 CM)	CX. 5
72382-08	LEUKOMED T PLUS 8 CM X 15 CM (PAD 4 CM X 11 CM)	CX. 5
72382-15	LEUKOMED T PLUS 10 CM x 20CM (PAD 5 CM X 16 CM)	CX. 5
72382-09	LEUKOMED T PLUS 10 CM X 25 CM (PAD 5 CM X 20,5 CM)	CX. 5
72382-10	LEUKOMED T PLUS 10 CM x 30 CM (PAD 5 CM X 25 CM)	CX. 5
72382-11	LEUKOMED T PLUS 10 CM x 35 CM (PAD 5 CM X 30 CM)	CX. 5