

# WEEKLY MEDICATION TRAY BLP-3100



## MAIN SPECIFICATIONS

- Weekly medication tray with 7 removable modules;
- Available in green and gray color;
- 7 removable containers;
- Adjustable and removable inner spaces;
- High capacity: ampoules tablets, blister;
- Lids that open in both directions;
- Tags in the containers sides;
- Washable up to 60 ° C. Do not use abrasive products;
- Stackable;
- Dimensions: 320 x 200 x 40 mm

Code	Description	Selling unit
BLP-3100	WEEKLY MEDICATION TRAY BLP-3100	CX. 12