

MEDICATION CART FOR 32 TRAYS



MAIN SPECIFICATIONS

- Medication cart for transportation of 32 BLP-3100 medication trays (not included);
- Capacity for 32 BLP-3100 medication trays (not included);
- 4 wheels with 125 mm, 2 with brakes and 2 simple (TENTE);
- Mini folding doors with lock (master key);
- Weight: 47,7 kg;
- ABS cover;
- Phenolic side;
- Inner guides of ABS / PV;
- Structure and base in galvanized steel epoxy;
- Accessories: rubbish (not included);
- Dimensions: Width: 717 mm; Length: 527 mm; Height: 1150 mm.

Code	Description	Selling unit
C-330A	MEDICATION CART FOR 32 TRAYS	CX. 1