

MEDICATION CART FOR 60 TRAYS



MAIN SPECIFICATIONS

- Medication cart for transportation of 60 BLP-3100 trays (not included);
- Capacity: 60 BLP-4100 trays (not included);
- Adapted walls/guides;
- 4 TENTE wheels with 125mm, 2 with brakes;
- 2 Doors/shutters with side opening (you don't to get down);
- Lock with key (digital lock is optional);
- Cover ABS / PC;
- Phenolic walls;
- ABS inner guides;
- Anodized aluminum profiles;
- Epoxy painted galvanized steel base;
- Dimensions: Width: 1116 mm; Length: 527 mm; Height: 1150 mm.

Code	Description	Selling unit
C-364	MEDICATION CART FOR 60 TRAYS TRANSPORTATION	CX.1