

# CART FOR THE TRANSPORTATION OF 60 TRAYS



## MAIN SPECIFICATIONS

- Medication cart for transportation of 60 BLP-4100 trays (not included);
- Capacity: 60 BLP-4100 trays (not included);
- Walls / adapted guides;
- 4 TENTE wheels with 125mm, 2 with brakes;
- 2 doors / shutters with side opening (do not have to get down), key lock (digital lock is optional);
- Cover ABS / PC;
- Phenolic walls;
- Interior ABS guides;
- Structure and base on galvanizado coated steel epoxy;
- Dimensions: Width: 1116 mm; Length: 527 mm; Height: 1410 mm

Code	Description	Selling unit
C-460	CART FOR 60 TRAYS TRANSPORTATION	CX. 1