

# MEDICATION CART FOR TRANSPORTATION OF 32 BLP-5100 TRAYS



## MAIN SPECIFICATIONS

- Medication cart for the transportation of 32 BLP-5100 trays (not included);
- Capacity: 32 BLP-5100 trays (not included);
- 4 TENTE wheels with 125mm, 2 with brakes;
- Door/ shutter with side opening (do not have to get down), lock with key (digital lock is optional);
- Cover ABS / PC;
- Phenolic walls;
- ABS inner guides;
- Anodized aluminum profiles;
- Epoxy painted galvanized steel base;
- Dimensions: Width: 717 mm; Length: 527 mm; Height: 1270 mm.

Code	Description	Selling unit
C-532	MEDICATION CART FOR 32 TRAYS TRANSPORTATION	CX. 1