

MEDICATION CART FOR TRANSPORTATION OF 36 BLP-6100 TRAYS



MAIN SPECIFICATIONS

- Medication cart for transportation of 36 BLP-6100 trays (not included);
- Capacity: 36 BLP-6100 trays (not included);
- Walls / adapted guides;
- 4 Wheels TENTE 125 mm, 2 with brakes;
- Side opening doors / shutters ;
- Locker with key (digital lock is optional);
- Cover: ABS / PC;
- Phenolic walls;
- ABS inner guides;
- Anodized aluminum profiles;
- Epoxy painted galvanized steel base;
- Dimensions: Width: 1110 mm; Length: 525 mm; Height: 1470 mm.

Code

Description

Selling unit

FACTORMED
Produtos Médicos & Hospitalares