

MEDICATION CART FOR THE TRANSPORTATION OF 16 BLP-6100 TRAYS



MAIN SPECIFICATIONS

- Medication car to transport 16 BLP-6100 trays (not included);
- Capacity: 16 BLP-6100 trays (not included);
- Walls / adapted guides;
- 4 Wheels TRY 125mm, 2 with brakes;
- 1 Door / shutter with side opening (do not have to download);
- Lock key (digital lock is optional);
- Cover in ABS / PC;
- Phenolic walls;
- ABS inner guides;
- Anodized aluminum profiles;
- Base in painted galvanized steel epoxy;
- Dimensions: Width: 720 mm; Length: 525 mm; Height: 1270 mm.

Code

Description

Selling unit

