

URINE BAG NON STERILE



MAIN SPECIFICATIONS

- Drainage bags for a very practical, hygienic and resistant urine collection.
- Composition: high quality PVC, free of heavy metals.
- Welds made by high frequency for maximum security.
- They have the top two reinforced holes to be able to hang the bag.
- Graduation, to measure the volume of urine.
- Universal connectors, adaptable to any type of civet.
- Protective cap on the connector.
- Connecting pipe resistant to disallow bottlenecks.
- Models available with or without output device in order to be emptied or not and with or without safety valve to control reflux or infections.

Code	Description	Selling unit
11011013	URINE BAG 1,5 L NON STERILE WITH NON RETURN VALVE AND 90 CM TUBE	CX. 1
11011016	URINE BAG 1,5 L NON STERILE WITH 90 CM TUBE	CX. 1
11011000	URINE BAG 2 L NON STERILE WITH NON RETURN VALVE, TAP AND 90 CM TUBE	CX. 1

