

## JELONET COTTON DRESSING



### MAIN SPECIFICATIONS

- Composed of knitted cotton fabric impregnated with paraffin.
- Soothes and protects the wound.
- It allows free passage of exudate.
- It can be used in combination with a secondary dressing and antiseptics.
- Suitable for: burns, leg ulcers, wounds with loss of skin, lacerations and abrasions.
- They are provided sterile in a variety of sizes.
- Soothes and protects the wound and allows free passage of viscous exudate.
- It has a soothing low tack which allows the wound to drain freely into the absorbent secondary dressing.
- Retains shape without fraying, open woven gauze dressings made have connected wires to be cut not fray.
- It can be combined with topical medication of choice.
- It does not have medicine and is ideal for topical use with antiseptic or antibiotic of choice.
- Protects injury preventing adhesion of the dressing shield and preventing the dehydration of the granulation tissue.
- It is indicated for the treatment of skin lesions clean, protecting and preventing rupture of the granulation tissue on dressing change. Suitable for frequent changes.

Code	Description	Selling unit
01003014	JELONET 5X5 CM	CX. 50
01003013	JELONET 10X10 CM	CX.10