

MULTIFUNCTIONS CART FOR TRAYS TRANSPORTATION



MAIN SPECIFICATIONS

- Multifunctional car to transport ISO 600x400 trays (not included);
- Robust and lightweight car, resistant to all chemicals used in a hospital setting;
- Walls / guides adapted to most of these trays in the market;
- 4 TENTE wheels with 125mm, 2 with brakes;
- Door / window shade with side opening (do not have to download);
- Lock key (digital lock is optional);
- Accessories and impeller adjustable in height;
- Higher plane of ABS / PC with flanges on three sides;
- Phenolic walls;
- ABS inner guides;
- Anodized aluminum profiles;
- Base in painted galvanized steel epoxy;
- Dimensions: Width: 840 mm; Length: 527 mm; Height: 1150 mm;
-

Weight: 36 kg.

Code	Description	Selling unit
CONVIMED-CMR1	MULTIFUNCTIONAL CART FOR TRAYS TRANSPORTATION	CX. 1