

UNIT DOSE MEDICATION CART FOR BIN-01



MAIN SPECIFICATIONS

- Multifunctional carts, ISO 600x400 for 42 BIN-01;
- 4 Wheels TENTE 125mm, 2 with locking;
- Shutter, key lock (digital lock optional);
- ABS/PC worktop;
- Phenolic walls;
- ABS support system;
- PVC shutter;
- Anodized aluminium profile;
- Epoxi painted galvanized steel base;
- Sizes: Width: 840 mm; Depth: 527 mm; Height: 1050 mm
- Weight: 34-36 kg

Code	Description	Selling unit
CONVIMED-CMU1	UNIDOSIS MEDICATION CART BIN-01	CX. 1