

UNIT DOSE MEDICATION CART FOR BIN-01/BIN-02



MAIN SPECIFICATIONS

- Robust car and light, resistant to all chemicals existent at hospitals;
- For 30 patients;
- Walls / guides adapted to most trays existing in the market;
- 4 TENTE wheels 125mm, two with brakes;
- Door/ shutter with side opening (do not have to downloadget down);
- Lock key (digital lock is optional);
- Accessories and impeller adjustable in height;
- Higher plane of ABS / PC with flanges on three sides;
- Phenolic walls;
- ABS inner guides;
- Anodized aluminum profiles;
- Base in painted galvanized steel epoxy;
- Dimensions: Width: 840 mm; Length: 527 mm; Height: 1150 mm;
-

Weight: 34-36 kg.

Code	Description	Selling unit
CONVIMED-CMU2	UNIDOSIS MEDICATION CART BIN-01/ BIN-02	CX. 1