

# UNIT DOSE MEDICATION CART FOR BIN-02



## MAIN SPECIFICATIONS

- Robust and lightweight car, resistant to all chemicals existent at hospitals;
- Walls / guides adapted to most of the trays in the market;
- 4 wheels 125mm TENTE, 2 with brakes;
- Door/ shutter with side opening (do not have to get down);
- Lock key (digital lock is optional);
- Accessories and impeller adjustable in height;
- Higher plane of ABS / PC with flanges on three sides;
- Phenolic walls;
- ABS inner guides;
- Anodized aluminum profiles;
- Base in painted galvanized steel epoxy;
- Dimensions: Width: 840 mm; Length: 527 mm; Height: 1150 mm;
- Weight: 34-36 kg.

Code	Description	Selling unit
CONVIMED-CMU3	UNIDOSIS MEDICATION CART BIN-02	CX. 1

FACTORMED  
Produtos Médicos & Hospitalares