

MULTIFUNCTIONAL CART WITH 5 DRAWERS



MAIN SPECIFICATIONS

- Robust car and light, resistant to all chemicals used at hospitals;
- 5 drawers: 4 drawers 600x400x100 mm; 1 drawer 600x400x200 mm;
- ISO modules not included;
- Front compact 13 mm thickness;
- Available colors: blue, red and white;
- Load capacity: 30 kg;
- Inner tray removable easily;
- Lock option with key in the drawers;
- 125 mm wheels, two with brake;
- Higher plane of ABS with flanges on three sides;
- Phenolic walls;
- ABS inner guides;
- Anodized aluminum profiles;
-

Structure epoxy steel;

- Dimensions: Width: 710 mm; Length: 525 mm; Height: 1050 mm;
- Weight: 36 kg.

Code	Description	Selling unit
D6040-5 0410	MEDICAL CART WITH 5 DRAWERS	CX. 1