

# MULTIFUNCTIONAL CART WITH 5 DRAWERS



## MAIN SPECIFICATIONS

- Robust and lightweight car, resistant to all chemicals existing at hospitals;
- 5 drawers: 1 drawer 600x400x50 mm; 3 drawers 600x400x100 mm; 1 drawer 600x400x200 mm;
- Available colors: blue, red and white;
- Capacity: 30 kg;
- Optional key that locks all drawers;
- ISO modules not included;
- Front in a compact 13 mm thickness;
- Removable inner tray easily;
- 125 mm wheels, two with brakes;
- Higher plane of ABS with flanges on three sides;
- Phenolic walls;
- ABS inner guides;
- Anodized aluminum profiles;
-

Structure in epoxy steel.

- Dimensions: Width: 710 mm; Length: 525 mm; Height: 1050 mm;
- Weight: 36 kg.

Code	Description	Selling unit
D6040-5 1310	MEDICAL CART WITH 5 DRAWERS	CX. 1