

## MEDICAL CART WITH 6 DRAWERS



### MAIN SPECIFICATIONS

- Robust and lightweight car, resistant to all chemicals existing at hospitals;
- 6 drawers: 1 drawer 600x400x50 mm; 5 drawers 600x400x100 mm;
- Front in a compact 13 mm thickness;
- 3 colors available: blue, red and white;
- Capacity / drawer: 30 kg;
- Optional key that locks all drawers stand;
- Possibility of inner tray extraction easily;
- 125 mm wheels, two with brakes;
- Higher plane of ABS with flanges on three sides;
- Phenolic walls;
- ABS inner guides;
- Anodized aluminum profiles;
- Structure in painted steel epoxy;
- Dimensions: Width: 710 mm; Length: 525 mm; Height: 1050 mm;
- Weight: 36 kg.

Code	Description	Selling unit
D6040-6 1500	MEDICAL CART WITH 6 DRAWERS	CX. 1