

ANESTHESIA CART



MAIN SPECIFICATIONS

- Robust and lightweight car, resistant to all chemicals existing at hospitals;
- 5 drawers: 1 drawer 600x400x50 mm; 3 drawers 600x400x100 mm; 1 of 600x400x200 mm drawer;
- ISO modules not included;
- Front in a compact 13 mm thickness;
- Available colors: blue, red and white;
- Capacity / drawer: 30 kg;
- Aluminum corners that allow the installation of accessories without prior hole;
- Possibility of inner tray extraction easily;
- 125 mm wheels, two with brakes;
- Higher plane of ABS / PC with flanges on three sides;
- Phenolic walls;
- ABS inner guides;
- Profiles in anodized aluminum;
- Structure in epoxy galvanized steel;
- Dimensions: Width: 710 mm; Length: 525 mm; Height: 1050 mm;

- Weight: 36 kg.

Code	Description	Selling unit
D-ANEST	ANESTHESIA CART	CX. 1