

TREATMENT CART



MAIN SPECIFICATIONS

- Treatment Cart, ISO 600x400. (ISO baskets not included);
- Upper removable ABS top cover with extended up rims on three sides;
- ISO drawers;
- Total opening 110%, telescopic ball bearing slides;
- Safe self-closing system, drawers;
- 4 Wheels TENTE 125mm 2 with locking;
- Optional 2 anti-static wheels;
- Ergonomic pushing handle;
- ABS/PC worktop;
- Phenolic walls;
- ABS support system;
- PVC shutter;
- Anodized aluminium profile;
- Epoxi painted galvanized steel base;
- Sizes: Width: 710 mm; Depth: 525 mm; Height: 1072 mm;
- Weight: 36 Kg.

Code	Description	Selling unit
D-TREAT	TREATMENT CART	CX. 1