

MULTIFUNCTIONS CAR FOR TRAYS TRANSPORTATION



MAIN SPECIFICATIONS

- Trays to transport ISO cabinets 600x400 (not included);
- Walls / guides adapted to most market trays;
- 4 Wheels TENTE 125mm, two with brakes;
- Shutter door with side opening (do not have to get down);
- Lock key (optional digital lock);
- Cover ABS / PC;
- Phenolic walls;
- ABS inner guides;
- Anodized aluminum profiles;
- Base in painted galvanized steel epoxy;
- Dimensions: Width: 840 mm; Length: 527 mm; Height: 1500 mm.

Code	Description	Selling unit
C-6040-1500	MULTIFUNCTIONAL CART FOR ISO 600 X 400 MM TRAYS	CX. 1

TRANSPORTATION

FACTORMED
Produtos Médicos & Hospitalares