

NURSING OPEN CART



MAIN SPECIFICATIONS

- Open carts with a flexible racking system designed for both angled and horizontal insertion of supports;
- Metallic structure;
- ABS Guides easy to move;
- Compact top;
- Wheels with 7,5 mm, 2 with brakes;
- The supports can be adjusted with 64 mm spaces;
- Changes and additions to the system are easy to make;
- Supports are designed with a stop-function - in the front and rear;
- The freestanding frames are made of 25x25 mm metal tubes and powder coated with epoxy;
- All LaPastilla modules fit into the LaPastilla open carts- adaptable to trays 600 x 400 mm with different heights (trays are not included);
- Dimensions: Width: 445 mm; Depth: 685 mm; Height: 1072 mm.

Code	Description	Selling unit
EM-04	NURSING OPEN CART	CX. 1