

MATTRESS LATEX WITH 12 CM ANT-BEDSORE AFRODITE DENSITY



MAIN SPECIFICATIONS

- Hospital mattresses 100% latex. The latex is a natural product very porous with optimal ventilation; hypoallergenic product that naturally prevents the growth of bacteria, fungi and mites inside;
- Firm and elastic material that allows a correct support of the spinal column and facilitates blood circulation;
- Articulated, reversible mattress with two equal sides;
- 12cm latex core, three densities available, depending on body weight or user preferences;
- It has a waterproof double elastic cover, liquids and gas-permeable, breathable allowing the breath of the interior, facilitating the hygiene of the mattress;
- Fire-resistant in accordance with BS 7175 1989, Crib 5;
- You can choose between three densities depending on body weight or user preferences:
 - **Normal:** Latex with density 60kg/ m³ for users weighing up to 70kg or preference for low consistency mattresses;
 - **Firm:** Latex with a density of 70kg/ m³, for users who weigh up to 85kg or preference for medium consistency mattresses;
 - **Extra Firm:** Latex with a density of 80kg/ m³, for users with body weight above 85kg or preference for high consistency mattresses.

- Dimensions: 190x90x12cm.
- Note: we accept orders for tailored mattresses for measures (by budget).
- Mattress Eros available- Dimensions: 190x90x15cm.

| Code | Description | Selling unit |
|---------------|---|--------------|
| 30009001 | MATTRESS LATEX WITH 12 CM ANT- BEDSORE AFRODITE DENSITY NORMAL UP TO 80 KG | CX. 1 |
| 30009002 | MATTRESS LATEX WITH 12 CM ANT- BEDSORE AFRODITE DENSITY FIRM + 80 KG | CX. 1 |
| 61 CL HC ER N | MATTRESS EROS 190X90X15 CM | CX. 1 |