

PAD ANTI- BEDSORE FLOATING NET



MAIN SPECIFICATIONS

- Sigel seat pad;
- Provides a larger contact area, which reduces the pressure per surface unit, one of the factors responsible for bedsores appearance;
- The Sigel liquid offers an unique sensation of coolness and confort;
- Made in an anti- slip base and comfortable upper material, waterproof;
- Made in an anti- slip base and comfortable upper material, waterproof;
 - Gel pad anti-bedsores for users up to 80 kg: 41x41x3,5cm;
 - Gel pad anti-bedsores horseshoe shape for users up to 80 kg: 41x41x3,5cm;
 - Gel pad anti- bedsores for users up to 100 kg: 42x42x4,5cm;
 - Gel pad anti-bedsores horseshoe shape for users up to 100 kg: 42x41x4,5cm;
 - Gel pad anti- bedsores for users up to 120 kg: 42x42x6cm.

Code	Description	Selling unit
30009003	PAD ANTI- BEDSORE FLOATING NET UP TO 80 KG	CX. 1
30009004	PAD ANTI- BEDSORE FLOATING NET UP +100 KG	CX. 1
30009011	PAD ANTI- BEDSORE FLOATING NET UP TO 100 KG	CX. 1

FACTORMED
Produtos Médicos & Hospitalares