

# ANTI BEDSORES MATTRESS W/ COMPRESSOR W/ REGULATION DOMUS I



## MAIN SPECIFICATIONS

- The anti-bedsores mattresses are known to be more effective preventing and complementing bedsores treatment;
- They are designed to reduce the risk of bedsores, even in persons who have a large skin sensitivity; its greater focus is the pressure relief and patient comfort;
- The new mattress is designed to be used as pillow-top mattresses pressure switches to prevent bedsores. This system is used to the treatment of grade I bedsores and supports a maximum weight of 100kg;
- The compressor is provided with a pressure regulator to suit the specific situation of each patient;
- With an almost imperceptible sound and its attractive design, we are truly talking about a choice that guarantees quality and an efficient treatment;
- Compressor:
  - Voltage: AC 230V 50Hz;
  - Air flow: 4 liters per minute;
  - Pressure range: 50 - 105 mmHg;
  - Cycle time: 9.6 minutes / 50Hz;

- Dimensions (W x L x H): 25 x 12.5 x 19,5cm;
- Weight: 1.9kg;
- Mattress:
  - Structure: 6.35cm with incorp extensions .;
  - Material (covers, cells and base): PVC with flame-retardant material;
  - Dimensions (W x L x H): 200 x 86 x 6,40cm;
  - Patient Weight: 100kg;
- Patches kit with glue is available.

Code	Description	Selling unit
30009009	ANTI BEDSORES MATTRESS W/ COMPRESSOR W/ REGULATION DOMUS I	CX. 1
30009015	PATCHES KIT WITH GLUE	CX. 1