

CIRCULAR CLOSED SEAT ANTI BEDSORES POLYURETHANE



MAIN SPECIFICATIONS

- The anti-bedsores seats are an asset in the dispersion of the loads that people are continually subject to and, sometimes, are areas of sores;
- Comfort tool;
- Material: polyurethane; blue fabric;
- Outside: bielastic cover, fireproof, impervious to liquids; It is breathable;
- Interior: polyurethane foam with a density of 50kg / m³;
- Outside diameter: 37cm;
- Inside diameter: 12 cm;
- Thickness: 7 cm.

Code	Description	Selling unit
30009012	CIRCULAR CLOSED SEAT ANTI BEDSORES POLYURETHANE	CX. 1