

# STATIC ANTI BEDSORES MATTRESS

## ATHENA 190x90 CM



### MAIN SPECIFICATIONS

- Mattress designed to minimize the pressure on the skin, in order to contribute to the prevention of decubitus sores;
- The top mattresses must be installed on a mattress in good conditions and with a smooth surface;
- Allows you to adjust the length of the mattress on which it is placed;
- Filling in buds hollow siliconized fibers, 100% polyester, not resinous fiber and coating 100% cotton fabric;
- Outer layer: bielastic, impermeable to liquids and permeable to gases, transpirável- allows breathing of the interior, facilitating the hygiene of the mattress. Fire-resistant according to BS 7175 1989, Crib 5. The tissue can be washed to a temperature of 95 °C with neutral detergent;
- Dimensions: 210x90x15cm.

Code	Description	Selling unit
30009013	STATIC ANTI BEDSORES MATTRESS ATHENA 190x90 CM	CX. 1